


# PHONOSORB® Desiccant Grades

PHONOSORB® molecular sieve beads are available in different particle sizes for different spacer bar systems and filling techniques:



Available in 25kg boxes or 150kg drums

- 0.5-0.9 mm / 35-20 mesh beads and automatic filling systems where the desiccant is pneumatically loaded into the profile through small injectors
- 0.8-1.6 mm / 20-12 mesh beads have been engineered for use in automatic spacer bar bending machines
- 0.8-2.0 mm/ 18-10 mesh beads with an optimum particle size for conventional filling of spacer bars

Distributed by:

**National  
Seal Systems  
Limited**

*Unit 210, Holly Road, Western Industrial Estate,  
Naas Road, Dublin 12  
Ph: 00353-1-4565170 Fax: 00353-1-4502371  
Email: [info@nationalseal.ie](mailto:info@nationalseal.ie) web: [www.nationalseal.ie](http://www.nationalseal.ie)*

Contact us for all your IG requirements